This Owner's Manual is provided and hosted by [Appliance Factory Parts](#).

# PANASONIC EP790 Owner's Manual

[Shop genuine replacement parts for PANASONIC  
EP790](#)

[Find Your PANASONIC Massage Lounger Parts - Select From 25 Models](#)

----- Manual continues below -----

# Operating Instructions

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## Manual de utilización

Household Massage Lounger

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
Model No. **EP790**

Modèle No. **EP790**

Número de modelo. **EP790**

**Panasonic**<sup>®</sup>

Before operating this unit, please read these instructions completely, and save for future reference.  
Il est recommandé de lire attentivement ce manuel avant d'utiliser l'appareil. Conserver ce manuel.  
Lea completamente estas instrucciones antes de utilizar esta unidad, y guárdelas para consultarlas en el futuro.

## TABLE OF CONTENTS

IMPORTANT SAFETY INSTRUCTIONS .....	3	MANUAL MASSAGE .....	13-15
SAFETY PRECAUTIONS .....	4-5	RECLINING THE CHAIR .....	16
GROUNDING INSTRUCTIONS .....	6	HOW TO EXTEND THE LEGREST .....	16
PARTS IDENTIFICATION .....	7	HEADREST .....	17
SPECIFICATIONS .....	7	HOW TO CLEAN .....	17
TURNING ON THE POWER .....	8	PLACE OF USE .....	17
CONTROLLER .....	9	MASSAGING HEADS RANGE MOVEMENT .....	18
TO START/STOP MASSAGE ACTION .....	10	TROUBLESHOOTING GUIDE .....	19
PROGRAMMED MASSAGE .....	11-12		

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### - Customer's Record -

Please note both the model and serial numbers of this product in the spaces provided. Retain this book, along with a copy of your purchase receipt, for future reference.

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Date Purchased

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# THANK YOU MESSAGE

Thank you for purchasing a Panasonic Massage Lounger.  
For optimum performance and safety, please read these instructions carefully.

## IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:

Read all instructions before using the massage lounger.

### **DANGER** - To reduce the risk of electric shock:

1. Always unplug this appliance from the electrical outlet immediately after using and before cleaning.

### **WARNING** - To reduce the risk of burns, fire, electric shock or injury to persons:

1. An appliance should never be left unattended when plugged in. Unplug from outlet when not in use.
2. Close and continuous supervision is necessary when this appliance is used by invalids or disabled persons.
3. Use this appliance only for its intended use as described in this manual.
4. Do not use attachments not recommended by the manufacturer.
5. This appliance should not be used by children. Keep children away from massage lounger. To avoid possible injury always return chair to an upright position and close the legrest carefully. Make sure nothing can be caught in the legrest mechanism.
6. Never operate this appliance if it has a damaged cord or plug, if it is not working properly or if it has been damaged. Call your nearest authorized servicer for repair.
7. Keep the cords away from heated surfaces.
8. Never drop or insert any object into any opening.
9. Do not use outdoors.
10. Do not operate where aerosol spray products are being used or where oxygen is being administered.
11. To disconnect, turn all controls to the off position, then remove plug from outlet.
12. Connect this appliance to a properly grounded outlet only.  
See Grounding Instructions.
13. To avoid the risk of electric shock or personal injury, do not attempt to open or disassemble any part of this massage lounger. Never insert hands in areas of moving parts. There are no user serviceable parts.

## SAVE THESE INSTRUCTIONS

# SAFETY PRECAUTIONS

## WARNING: To avoid risk of injury

1. This product should not be used by pregnant women, persons who are ill or persons with back, neck, shoulder or hip pain, or conditions of any kind, without consulting a physician prior to use. Do not use on any area that has received medical treatment without consulting a physician prior to use.
2. Never put any part of the body between the massage heads. The squeezing action of the massage heads may cause injury. Never put any part of the body between the arm rest and the chair back.  
Do not place hands in mechanical portion of legrest located under the seat.
3. Do not use without clothing.  
Thin clothing increases the effectiveness, but exposing the skin directly to the massager can irritate the skin.
4. Do not use for massaging the head, abdomen, elbows or knees.
5. Do not stretch the back muscles for more than approx. 15 minutes at one time.  
To avoid overstretching of muscles and resulting discomfort when first using the chair, do not use for more than 5 minutes. You can gradually increase duration and intensity as you adjust to the chair.  
Do not massage any one point for longer than 5 minutes at a time.  
Excessive massaging can overstimulate the muscles and nerves and result in an adverse effect.  
To prevent excessive massaging, do not use the massager for more than a total of 15 minutes per session. After approx. 15 minutes, the massage heads will return to the "restore" position and stop.
6. While using the massager, if you start feeling sick, or if the massage seems painful, stop use immediately.  
While using the massager, if the pressure seems too strong, the movement of the massage heads can be stopped by pressing the stop button.
7. When the chair back is in a reclined position, do not sit on it, or sit with your legs resting on the headrest.  
Do not sit on the legrest. The chair may fall over.
8. Do not allow children or pets to play on, or around, the chair, especially during operation.
9. Do not remove the back cushion without first reading instructions on page 17.

## **CAUTION: To avoid damage to the chair**

1. Do not sit on, or drop, the controller. Do not pull the controller cord or the power cord. Be careful not to catch the cord of the controller between the arm rest and the chair back.
2. Do not treat the chair roughly, such as turning the chair on its side, turning the chair over and standing on it.
3. Do not treat the covering roughly.
  - Keep sharp or pointed objects away from the cloth of the chair. Be careful not to drop lit ashes or matches on the chair.
  - Exposure to direct sunlight can cause fading or a change in colour of the covering.
4. If the massager functions abnormally, immediately turn off the power and have the massager checked by an authorized servicenter.
5. Do not attempt to repair the massager yourself.
6. Do not repeatedly operate any switch at short intervals. It may cause the switch to malfunction.
7. Some sounds may be heard coming from the massager while in use. This is due to the structure of the massager. This is a normal condition.
8. Make sure that there are no obstacles behind the chair when reclining it. The chair can be reclined to a maximum angle of 170°. If the back hits a wall or pillar, it may malfunction. Be sure to have ample space behind the chair.
9. When excessive pressure is applied to the massage heads, they may stop moving for safety. When moving the heads up from the lowered position, the body weight may stop the heads movement. In such case, slightly raise the body and allow the heads to move up.
10. After each massage, slide the power switch, which is located on the side of the chair, to “off” position and turn its lock switch to “lock” position.  
**NOTE:** Maximum user weight is 264 lbs (120 kg).

## **CARING FOR THE MASSAGER**

- Always unplug chair before cleaning.
- If the covering becomes soiled, dampen with liquid soap and brush the surface. Wipe with a cloth dampened with water. Allow to dry naturally. If the plastic portions become soiled use lukewarm water, with liquid soap on a soft cloth to clean the surface.
- Do not use thinner, benzene or other solvent type cleaners.

# GROUNDING INSTRUCTIONS

This product must be grounded. If it should malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock.


This product is equipped with a cord having an equipment-grounding conductor and a grounding plug. The plug must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances.

## DANGER

Improper connection of the equipment-grounding conductor can result in a risk of electric shock. Check with a qualified electrician if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product if it will not fit the outlet. Have a proper outlet installed by a qualified electrician.

This product is for use on a nominal 120-volt circuit and has a grounding plug that looks like the plug illustrated in Figure A. Make sure that the product is connected to an outlet having the same configuration as the plug. No adaptor should be used with this product.

**Figure A**

## PARTS IDENTIFICATION

### Massage lounger

#### Back cushion

- Attaches to the chair with a zipper.
- Remove if the massage feels weak.

#### Headrest

- Attaches to the back cushion with a zipper.

#### Massage heads

(Underneath chair cushion)

#### Legrest


- With leg massager (Leg massager can be used using the legrest extender lever.)

#### Armrest

#### Controller storage pocket

(Controller can be stored in the controller pocket (left side).)

## SPECIFICATIONS


Power source	: 120 V AC, 60 Hz
Power consumption	: 230 W
Kneading speed	: Approx. 28 times/min.
Tapping speed	: Approx. 500 times/min. (per side)
Rolling massage speed	: Approx. 1 cycle every 33 sec.
Massage width	: Shoulder/lower back section : Approx. 2 15/16 in. (75 mm)
Back rolling width/tapping width	: Narrow : Approx. 3 5/16 in. (85 mm) Wide : Approx. 4 5/16 in. (110 mm)
Massage heads up/down travel	: Approx. 23 1/16 in. (600 mm)
Regional back rolling	: Automatic repetition within approx. 4 12/16 in. (120 mm) range
Intensity adjustment	: Adjusts massage head protrusion steplessly within approx. 1 3/4 in. (45 mm) range
Shoulder position adjustment	: 7 steps
Automatic shut-off	: Approx. 15 min.
Leg massage - vibration	: Approx. 3200 cycles/min. Shut-off timer : Approx. 30 min.
Dimensions (H x W x D)	: • Not reclined and legrest stored : 41 11/32 x 29 23/32 x 41 11/32 in. (1050 x 755 x 1050 mm) • Reclined and legrest extended : 24 13/32 x 29 23/32 x 68 57/64 in. (620 x 755 x 1750 mm)
Reclining angle	: Approx. 127° to 170°
Weight	: 86 lbs. (56 kg)
Accessories	: Back cushion, seat cover, headrest

## TURNING ON THE POWER

First read, and familiarize yourself with, the safety precautions.


- Plug the power cord into the inlet located at the rear of the chair.  
Plug the power cord into an AC outlet.  
Be careful not to pinch the cord under the leg frames of the chair or any other items.

1. Using a flat screwdriver or similar object, insert into the lock switch and turn to "open".
2. Move the power switch to "on".
  - When the operating lock switch is pointing toward "lock", the power switch cannot be moved to the "on" position.

### After each use.


- Be sure to turn the power switch to "off".
- To prevent children from using this unit, safety lock the power switch by moving the operating lock switch to the "lock" position.
- As a further caution, unplug the power cord from the outlet after each use.

## BEFORE SITTING DOWN


After checking the positions of the massage heads with your hands, sit down slowly in the center of the Massage Lounger.

- The massage heads should be spread apart near the top of the Massage Lounger backrest (home position).
- If the massage heads are not in the position described above, do not sit down in the Massage Lounger. Press the "start/reset" button on the control panel twice. Wait for the massage heads to reach the home position before sitting down.


# CONTROLLER

## PROGRAMMED MESSAGE

## MANUAL MESSAGE


## TO START/STOP MESSAGE ACTION

Before using, first read all the instructions including the warnings and safety precautions as well as the section on programmed and manual massage.


### <BEGINNING OPERATION>

1. To operate the massager, press the "start/reset" button.
  - Check to be sure the operation indicator lamp is lit.


2. If you wish to stop the massage immediately, press the "stop" button. The massage heads will stop where they are.

- Check to be sure the operation indicator goes out. To restart massage press "start/reset" button.


3. After final use, press the "start/reset" button.

- The massage heads will go up, spread out and stop in the restore position. With the heads in this position you can sit comfortably, using the massager as a regular chair.
- Be sure to store the controller in the pocket on the side of the left arm rest when not in use.

4. Do you want programmed or manual operation?

- If you wish to select a programmed massage, slide the cover down fully and press buttons.
- If you wish to select the manual massage, slide the cover upwards fully and press buttons.

### Automatic timer for massage lounger

- To prevent excessive massaging there is a built-in timer which, after approx. 15 minutes of operation, will automatically return the massage heads to the "restore" position and stop their motion.
- Turning on the operation switch while the massage heads are moving to the "restore" position will reset the time for another 15 minutes (approximately).
- If the programmed massage button is pushed, the unit will stop after finishing the program.


## PROGRAMMED MESSAGE

**First read, and familiarize yourself with, the safety precautions.**

### <HOW TO PRESET THE MESSAGE HEAD POSITION>

1. Press the "start/reset" button with cover slide down fully.
2. Set the message heads adjustment switch to position "7".  
Then adjust the message heads to your shoulders while moving the adjustment switch.  
The message heads operation lamp flashes while moving the message heads.
3. You can adjust the message heads position while massaging.

**NOTE:** When you find your favourite position, memorize the number which is printed near the switch.


### <HOW TO SELECT THE PROGRAMMED MESSAGE COURSE>


1. If you wish to massage your whole back, press (A).
  2. If you wish to massage your neck and shoulders mainly, press (B).
  3. If you wish to massage your lower back, press (C).
- Message will operate automatically for 15 minutes after pressing the button.

**NOTE:** Check to be sure the course indicator lamp is lit.


- If you wish to change the massage course, press any button you want.  
You may change the massage course at any time during operation without harm.  
Your new instructions will override any previous ones.

<whole back>

<neck shoulder>

<lower back>

### <GENTLE MODE FOR PEOPLE WHO FEEL MESSAGE IS TOO STRONG>

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- Gentle mode indicator lamp is lit.  
The entire message changes to a less intense, more gentle message.  
(Course contents do not change)
- Gentle mode can be selected even in the middle of a programmed message course.
- Pressing the gentle mode button again will cancel gentle mode.


**After the automatic course has finished, the message heads will automatically move to the "restore" position and stop.**

- To perform message again, press the "start/reset" button on the controller again.


# PROGRAMMED MASSAGE

## Details of Automatic Course


### Whole back course

### Neck and shoulder course


### Lower back course

The above diagrams have been abbreviated to enable easy understanding of the massage action.

## MANUAL MASSAGE

First read, and familiarize yourself with, the safety precautions.

1. Slide the cover upward fully.
2. Push the massage mode button.
  - Massage operation can be changed by pushing the massage mode button progressively as follows.

1. Shiatsu/upward massage
2. Shiatsu/downward massage
3. Tapping
4. Rolling/stretching massage/whole back
5. Rolling/stretching massage/regional
6. Tapping and rolling/stretching massage/whole back
7. Tapping and rolling/stretching massage/regional

*Try out the various modes and see which suits you best!*

You can change the area by pressing and buttons.

### Shiatsu/upward/downward/massage

- ① Choose the “upward” or the “downward” mode and the massage begins.
- ② To adjust to the desired position, use the “up”(▲) or the “down”(▼) buttons.
 

**NOTE:** Do not massage any one point for more than 5 minutes. Do not stretch the back muscles for more than approx. 15 minutes at one time. A continuous massage in the same area may result in soreness.
- ③ Use the intensity adjustment buttons to adjust the intensity of the massage.

### Tapping massage

- ① Choose the “tapping” mode and the massage begins.
- ② Use the “narrow”(▶◀) or the “wide”(◀▶) adjustment button to adjust the width of the massage heads. To raise the massage heads, press the “up”(▲) button. To lower them, press the “down”(▼) button.
- ③ Use the intensity adjustment buttons to adjust the intensity of massage.

### Rolling/Stretching massage Whole back/regional

- ① Choose the “rolling” or the “regional rolling” mode.
  - The massage heads move to the narrowest position and start the rolling massage.
- ② Use the “narrow”(▶◀) or the “wide”(◀▶) select button to select the width of the massage heads.
- ③ To raise the massage heads, press the “up”(▲) button. To lower them, press the “down”(▼) button.
  - When you wish to stretch only a certain part of your back, choose the “regional rolling” mode .

All functions which include “rolling”, are available for the whole back (indicated by a red lamp) or a limited section of the back (indicated by a green lamp).

## MANUAL MASSAGE


### Tapping and rolling massage

- ① Choose the “tapping rolling” or the “regional tapping rolling” mode.
  - The massage heads move to the most gentle position and the narrowest position and start the tapping and rolling massage.
- ② Use the “narrow”(▶◀) or the “wide”(◀▶) adjustment button to adjust the width of the massage heads.
- ③ To raise the massage heads, press the “up”(▲) button. To lower them, press the “down”(▼) button.
  - When you wish to stretch only a certain part of the back, choose the “regional tapping rolling” mode. You can change the area by pressing(▼) and (▲) buttons.

#### Adjustability for Manual Massage

Manual Mode	Up/down Adjustment	Intensity Adjustment	Width Adjustment
upward/downward	○ 23 1/16" (600 mm)	○ 1 3/4" (45 mm)	×
tapping	○ 23 1/16" (600 mm)	○ 1 3/4" (45 mm)	○ 3 5/16" or 4 5/16" (85 or 110 mm)
rolling/tapping rolling (stroke: 600 mm)	○ 23 1/16" (600 mm)	○ 10/16" (15 mm)	○ 3 5/16" or 4 5/16" (85 or 110 mm)
regional rolling/regional tapping rolling (stroke: 600 mm)	○ 23 1/16" (600 mm)	○ 10/16" (15 mm)	○ 3 5/16" or 4 5/16" (85 or 110 mm)

### Massage intensity adjustment

- Hold your finger on the intensity adjustment button () to choose the massage intensity desired.
- When you adjust the massage intensity, the massage heads position may change. Adjust the massage heads position again.
- The massage heads stop automatically when they reach the maximum/minimum last intensity positions. Release your finger from the button and operation will resume.

## Message heads up-down and width adjustment ▲/▼ ▶◀/◀▶

### Adjust the vertical position of the message heads.

- Use the up or down adjustment button to adjust the up or down direction of the message heads.  
To raise the message heads, press the “up”(▲) button. To lower them, press the “down”(▼) button.
- The message heads will move while holding your finger on the button.
- When the message heads have reached the desired position, release the button and the heads will stay in that position.

### Select the width of the message heads.

Use the narrow (▶◀) or wide (◀▶) buttons to select the width between the message heads (available only in the tapping massage and rolling, tapping rolling).

**NOTE:** If the message heads are forcibly locked, they will stop for safety and the controller panel will warn you with blinking lamps.  
Raise your back slightly and press the “start/reset” button.

## RECLINING THE CHAIR


### First read, and familiarize yourself with, the safety precautions.

To avoid possible tipping of the chair, do not sit, or lean, on backrest while chair back is in a reclined position. The angle of the chair back can be changed so the body weight and the massage heads will work together to determine the strength of the massage.

Try changing the angle of the chair back to find a comfortable massage at a strength that is most suitable for your own comfort. Towels may be used in conjunction with the back cushion to lessen intensity, if desired.

### Leaning back

1


Press and hold 

While the button is held pressed, the legrest will extend and the backrest will recline.

**NOTE:** Position the chair at least 17 3/4 in. (45 cm) away from any wall to ensure adequate reclining space.

2

Release the button when the desired position has been reached.

### Returning back to upright after reclining

Press and hold 

- The legrest will be returned to its normal position (store) and the backrest will be raised.
- After use, always be sure to return the legrest and backrest to their normal positions (legrest stored and backrest raised).

## HOW TO EXTEND THE LEGREST


Grasp the lever.

**To extend: Tilt the lever back.**

- Tilt the lever back as far as it will go.

**To store: Lift the lever.**


- When extending or storing the legrest, lift your legs off of the legrest and operate the lever.

### NOTE:

- Keep children and pets away from the massage lounger, when lowering the legrest, to prevent injury.
- Although the legrest is moved automatically according to the movement of the Massage Lounger backrest during reclining, the legrest can also be raised and lowered independently using the side lever.
- The maximum height (angle) of the legrest when using powered reclining and when using the side lever may vary slightly.
- The legrest cannot be stored completely if the Massage Lounger backrest is not raised completely. (The same condition occurs for lever operation.)
- Be careful of the lever during powered reclining. It will move automatically.
- Since the legrest will be locked in position when it is raised using the side lever, it cannot be lowered using the reclining buttons. Always use the side lever to lower it if it has been raised using the side lever.

## HEADREST

- Fasten the back cushion to the chair and the headrest to the back cushion.
- When the massage function will be used, remove the headrest or flip the headrest over the back of the chair before starting the massage.
- If the massage seems weak, remove the back cushion.

## HOW TO CLEAN

- Wipe the cover with a cloth dampened with water and mild soap. Then, wipe with plain water and let dry.
- Do not use bleach or solvents such as benzene or alcohol. Chemicals can discolor or damage the material.

## PLACE OF USE

- Do not use the Massage Lounger where it can be exposed to direct sunlight for long periods of time. Over exposure to harsh light can discolour or harden the cover.
- Do not use the Massage Lounger in hot or humid places or places where it may become wet from rain, splashing or dripping water, etc. Water and heat can cause discoloring and mildew.


**NOTE:** This massage lounger is capable of providing a vigorous “Shiatsu” massage. The rolling and tapping are very strong, and that massage is not for people who have not previously experienced a shiatsu massage. For that reason, a foam cushion is provided to cover, and soften, the effect of the massage. The pad can be removed for a full shiatsu massage, but care should be taken. The first time you try it -- if you decide to -- sit in a relaxed, not tense, manner and do not try it for more than two minutes. After that you should be able to decide whether or not you prefer a massage with the pad in place.

## MASSAGING HEADS RANGE MOVEMENT

The width between the massage heads can be adjusted for tapping, rolling, regional rolling, tapping rolling massage and regional tapping rolling massage. The heads can be set at different widths for a wide range, giving you a massage any place you like from the neck to the lower back.


### <UPWARD, DOWNWARD MASSAGE>

Massage range width :  $2 \frac{15}{16}$  in. (75 mm) (width cannot be adjusted)


### <TAPPING, ROLLING, REGIONAL ROLLING, TAPPING ROLLING AND REGIONAL TAPPING ROLLING>

Massage range width :  $3 \frac{5}{16}$  in. or  $4 \frac{5}{16}$  in. (85 mm or 110 mm)

## TROUBLESHOOTING GUIDE

Before requesting service for this unit, check the table below for the possible cause of the problem. Some simple checks or minor adjustment on your part may eliminate the problem and restore proper operation. If you are in doubt about some of the check points, or if the remedies indicated in the chart do not solve the problem, refer to the directory of authorized servicers (see the servicer list included) to locate a convenient servicer, or consult your authorized Panasonic dealer for instructions.

Problem	Probable cause(s)	Suggested Remedy
Massage does not start.	The power cord is not plugged in.	Plug the power cord into an AC outlet (120 V).
	The power source switch is at the "off" position.	Slide the power switch to "on".
All lamps on controller flash and massage heads Stop.	If the massage heads are forcibly locked, they will stop for safety.	Raise your back slightly and push the start/reset button.
Kneading massage cannot be changed immediately to a rolling massage.	The massage heads must move to the narrowest position for a rolling massage.	Wait until the massage heads reach the narrowest position.
Massage heads do not reach your shoulders.	Your shoulders are too high.	Shift your waist forward and recline the chair back.
It seems like the massage heads are at different heights.	<p>This is because an alternating tapping mechanism is used. It is not a malfunction.</p>  <p>Approx. 13/32 in. (1 cm) at most</p> <p>Approx. 13/32 in. (1 cm) at most</p> <p>When stopped, the misalignment is approx. 1 cm at most.</p>	

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